

Taq-Purple DNA polymerase

(Catalogue number T107, T108, T109)

rev. 07/2025

Description

Taq-Purple DNA polymerase is an alternative to Taq DNA polymerase Unis (Cat. No. T037-T 039). In this product, concentration of Taq DNA polymerase is reduced to 1U/μl and storage buffer is supplemented with inert Dye for easy loading and control of the presence of the enzyme in PCR.

Inert red dye

The dye does not interfere with PCR and during electrophoresis in agarose gel it migrates faster than oligonucleotide primers. Therefore, the dye does not interfere with quantification of DNA amplicons.

Taq DNA polymerase

Taq DNA polymerase is a thermostable enzyme isolated from Thermus aquaticus. The enzyme catalyzes synthesis of complementary DNA strand in the 5'→3' direction and also possesses a 5'→3' exonuclease activity. During amplification of DNA fragments, Taq polymerase adds at 3' end an adenosine overhang. This can be utilized for cloning of PCR-generated DNA fragments. Advantage of the enzyme is its high processivity (amplification of 1000 base pairs takes < 1 min). Disadvantage is that it lacks a 3'→5' exonuclease proofreading activity and this accounts for high error rate [about 1 error to 10⁵ - 10⁶ base pairs (bps)]. The major usage of the enzyme is in diagnostic analysis for amplification of DNA fragments up to 5000 bps.

Technical data

Components and packaging

- Taq-Purple DNA polymerase is supplied at a concentration 1 U/μl. Basic packaging is 1 tube with 500 U/500 μl (**T107**), 5 tubes with 500 U/500 μl (**T108**) or 10 tubes with 500 U/500 μl (**T109**).
- Each tube of Taq DNA polymerase 1.1 is accompanied by a tube with 10x concentrated PCR Blue buffer containing MgCl₂ (1.5 ml).

If different concentration of MgCl₂ is required, a tube with 10x concentrated PCR Blue Buffer without MgCl₂ (1.5 ml) and a tube with 25 mM MgCl₂ (0.5 ml) should be ordered (Cat. No. **T059**).

Storage

- At temperature -20°C ± 5°C. Material can be repeatedly defrosted.

Composition

- Storage buffer: 20 mM Tris-HCl (pH 8.0 at 25°C), 100 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.5% Nonidet P-40, 0.5% Tween 20, inert red dye, stabilizers, 50% glycerol.
- 10x conc. PCR Blue Buffer: 750 mM Tris-HCl, pH 8.8 (at 25°C), 200 mM (NH₄)₂SO₄, 1% Tween 20, 25 mM MgCl₂.

Activity

- One unit is defined as the amount of enzyme, which catalyzes incorporation of 10 nmol dNTPs within 30 min at 72°C into trichloroacetic acid precipitable material. Reaction conditions are as follow: 10 mM Tris-HCl (pH 8.8 at 25°C), 50 mM KCl, 0.1% Triton X-100, 1.5 mM MgCl₂, 200 μM dATP, dCTP, dGTP and [α-³²P]dTTP, 50 μg/ml denatured cDNA, 0.5 μM primer and 0.2 – 0.5 U of enzyme in the volume of 50 μl.

Purity and quality control

- Purity of Taq-Purple DNA polymerase is tested by SDS-PAGE. Enzyme migrates as a major band of 94 kDa. Material is nuclease free.
- Each batch of Taq-Purple DNA polymerase is tested for its ability to amplify DNA fragment from mammalian genomic DNA by PCR. The results are verified by electrophoresis in agarose gel in the presence of ethidium bromide; only DNA band of the expected size is present.

Taq-Purple DNA polymerase

(Catalogue number T107, T108, T109)

Protocol

Basic protocol

The protocol described below can be used for routine PCR. However, in some cases reaction conditions must be optimized, mainly annealing temperature and concentration of MgCl₂.

1. In thin-wall test tubes the following components are mixed¹:

	PCR in 50 µl	Final concentration
10x PCR Blue Buffer with 25 mM MgCl ₂ ²	5 µl	1x react. buffer with 1.5 mM MgCl ₂
PCR dNTP mix (10 mM each) (Cat. No. P041)	1 µl	0.2 mM dNTP each
5' primer (50 µM)	0.5 µl	0.5 µM
3' primer (50 µM)	0.5 µl	0.5 µM
Taq-Purple DNA polymerase (1U/µl)	2.5 µl	2.5 U (0.05 U/µl)
Template DNA (1 ng/µl - 1 µg/µl)	1 µl	0.02 ng/µl – 0.02 µg/µl
PCR H ₂ O (Cat. No. P042)	39.5 µl	

¹ When more DNA samples are tested with the same primers, it is convenient to prepare so called Master Mix in which the volumes of individual components are multiplied by the number of DNA samples tested. Then, 49 µl aliquots of the Master Mix are distributed into the PCR tubes, followed by addition of tested DNA (1 µl) into each tube.

² If unsatisfactory results are obtained, PCRs with different concentrations of MgCl₂ should be performed (see below).

2. Samples are homogenized and spun down.

3. If cycler without heating lid is used, PCR oil 25 µl (Cat. No. P043) is added.

4. Thermal cycler is programmed according to the manufacturer's instructions. A typical cycling program is:

	Temperature	Time	Number of cycles
Initial denaturation	94°C	3 min	1
Denaturation	94°C	30 s	25-35
Annealing of primers	55-68°C ¹	30 s	
Extension	72°C	1 min per 1 kb	
Final extension	72°C	10 min	1
Cooling	4°C		

¹ Should be determined experimentally; usually 5°C below melting temperature of the primers.

5. After PCR, samples are mixed with loading buffer (Cat. No. P048, P066 or P065) and analyzed by electrophoresis in agarose gel in the presence of ethidium bromide. Alternatively, the samples can be stored at -20°C.

Optimization of MgCl₂ concentration

MgCl₂ at a 2.5 mM final concentration is suitable for most PCRs with Taq-Purple DNA polymerase. However, if amplification of nonspecific DNA fragments is observed, optimal Mg²⁺ concentration for given PCR should be determined. To this end 10x PCR Blue Buffer without MgCl₂ and 25 mM MgCl₂ can be ordered (Cat. No. T59).

1. Preparation of Master Mix **without MgCl₂** by mixing the following components:

10x reaction buffer without MgCl ₂	40 µl
PCR dNTP mix (10 mM each)	8 µl
5' primer (50 µM)	4 µl
3' primer (50 µM)	4 µl
Taq-Purple DNA polymerase (1U/µl)	20 µl
Template DNA (1 ng/µl - 1 µg/µl)	8 µl
PCR H ₂ O	236 µl
Total volume	320 µl

2. Master mix is thoroughly mixed, centrifuged briefly and 40 µl aliquots are distributed into 7 PCR tubes.

3. MgCl₂ and PCR H₂O is added into PCR Master mixes as follows:

Tube No.	25 mM MgCl ₂	PCR H ₂ O	Final MgCl ₂ concentration
1	2 µl	8 µl	1.0 mM
2	3 µl	7 µl	1.5 mM
3	4 µl	6 µl	2.0 mM
4	5 µl	5 µl	2.5 mM
5	6 µl	4 µl	3.0 mM
6	8 µl	2 µl	4.0 mM
7	10 µl	0 µl	5.0 mM

4. PCR is performed as above and the samples are analyzed by electrophoresis in agarose gel in the presence of ethidium bromide. Optimal concentration of MgCl₂ for given PCR is thus determined.

Cat. No.	Product name and specification	Amount
T107	Taq-Purple DNA polymerase	500U
T108	Taq-Purple DNA polymerase	5x 500U
T109	Taq-Purple DNA polymerase	10x 500U
T059	10x conc. PCR Blue Buffer without MgCl ₂ +MgCl ₂	1.5 ml + 0.5 ml

