

## PCR Ethanol, 96%

96% ETHANOL IN PCR H<sub>2</sub>O  
(Cat. No. P044)

rev. 01/2022

### Description

PCR Ethanol, 96% is 96% ethyl alcohol in ultrapure water (18 Mohm.cm, ultrafiltered), from which RNases were removed by the treatment with diethylpyrocarbonate followed by autoclaving.

### Technical data

#### Storage

- Store at temperature 4°C ± 3°C.

#### Packaging

- 100 ml in glass bottle with screw cap.

#### Quality control

- Each batch of PCR Ethanol, 96% is tested in RT-PCR.

Cat. No.	Product name and specification	Amount
P054	PCR Ethanol, 96%	100 ml

