

## CERTIFICATE OF ANALYSIS

**Product:** DNA marker 155-970 Coloured\_

**Catalog No:** C109, C111

**Lot No:** C109072028

**Date of Expiration:** 07/2028

**Concentration:** 1 µg DNA/12 µl buffer (25 mM potassium acetate, 10mM Tris acetate and 25 mM EDTA).

**Storage temperature:** -20 ± 5°C

**Test:** 10 µl of DNA marker 155-970 Coloured was directly loaded into wells of 1.2% agarose gel, which contained ethidium bromide (1 µg/ml) and 1x TBE or TA buffer. The bands under UV light should be visible (970 bp at the top).

**Result of the test:** passed

FOR RESEARCH USE ONLY

APPROVED DATE: 25.8.2025

Quality control: Hana Těšitelová